



Barrier contraceptives are making a strong return as enlightened couples seek hormone-free methods of birth control.

Increased momentum towards alternative lifestyles and healthier, natural living has brought this classic contraceptive option back into relevance.

The diaphragm forms a barrier between the sperm and the cervix and is used in conjunction with a contraceptive gel or spermicide.

Selecting the correct size is essential to maintaining contraceptive efficacy. A fitting kit containing the five most popular mid-range sizes is available to assist the healthcare professional in assessing a patient's correct size. The fitting kit can be autoclaved.

The size assessment should be repeated:

- Annually
- After pregnancy
- After weight loss/gain of more than 5kg

Our Milex brand of diaphragms is made in the USA. They are FDA approved and CE marked. They should be replaced two years after use, after any type of vaginal infection or if indication of wear is apparent.

Milex Diaphragm Key Points:

94% effective with perfect use.

Omniflex Wide-Seal rim can be folded in any direction.

Durable and latex-free.

Omniflex Diaphragm Sizing:

Milex Omniflex 55mm:
EAN: 5060244481985

Milex Omniflex 60mm:
EAN: 5060244480100

Milex Omniflex 65mm:
EAN: 5060244480117

Milex Omniflex 70mm:
EAN: 5060244480124

Milex Omniflex 75mm:
EAN: 5060244480131

Milex Omniflex 80mm:
EAN: 5060244480148

Milex Omniflex 85mm:
EAN: 5060244480155

Milex Omniflex 90mm:
EAN: 5060244480162

Milex Omniflex 95mm:
EAN: 5060244481992

Clinical Fitting Kit:
EAN: 5060244480544

Country of Origin:
United States

